

# RNase/DNase Decontaminant

CAT No. APL-1089

## DESCRIPTION

RNase/DNase Decontaminant is a non-alkaline, non-corrosive and non-carcinogenic cleaning solution which is highly active against nuclease contamination. In addition, RNase/DNase Decontaminant can effectively degrade Ethidium Bromide (EtBr) as quickly as 5 minutes in non-fluorescent and non-mutagenic environments. The solution contains a surfactant plus an agent that destroys RNase/DNase activity and degrades EtBr. RNase/DNase Decontaminant has been demonstrated to inactivate greater than 20 µg of RNase A dried powder onto the bottom of a microcentrifuge tube. RNase/DNase Decontaminant is stable and heat resistant. Storage at 65°C for 2 weeks will not reduce the quality of the product.

## CONTENT

	APL-1089
RNase/DNase Decontaminant	500 mL

## STORAGE

RNase/DNase Decontaminant should be stored at room temperature; at lower temperature, a precipitate may form which is easily brought into solution at 37°C.

## APPLICATION

RNase/DNase Decontaminant is ready-to-use for eliminating Nuclease and EtBr from any surface including the interior of microcentrifuge tubes. By following the simple decontamination protocol below, RNase and DNase is completely inactivated and removed.

## PROTOCOL

### To decontaminate pipette

Following manufacturers' instructions; remove the shaft from the pipette and then remove seals and gaskets from the shaft. Spray the shaft and stand for one minute in RNase/DNase Decontaminant, rinse the shaft thoroughly with water, let dry and assemble.

### To decontaminate laboratory surface

Spray RNase/DNase Decontaminant directly to the lab surface. Wipe thoroughly with a paper towel, rinse with water and dry with a clean paper towel.

**To decontaminate laboratory apparatus**

Generously apply RNase/DNase Decontaminant to a paper towel and wipe all exposed surfaces of the apparatus thoroughly. Rinse with water and dry with a clean paper towel. To clean small parts, briefly soak them in RNase/DNase Decontaminant, rinse with water and dry.

**To decontaminate plastic and glass vessel**

Spray RNase/DNase Decontaminant to enable coating the entire surface of the vessel or electrophoresis tank by swirling. Discard the solution and rinse vessels thoroughly two times with distilled water.

**NOTES:** *Prolonged contact of RNase/DNase Decontaminant with skin may cause irritation; therefore, gloves should always be worn.*

**PRODUCT USE LIMITATION**

Research use only.