
PRODUCT SPECIFICATION

Product Name: Agarose (Low EEO)

CAT Number: APL-0002

CAS Registry Number: 9012-36-6

Specification	
Appearance	White to off-white powder
DNase/ RNase/ Proteases	Not detectable
Electroendosmosis	0.09~0.13 (EEO)
Gel strength (1.5%)	$\geq 2300 \text{ g/cm}^2$
Gel point (1.5%)	$36 \pm 1.5^\circ\text{C}$
Melting Point	$88 \pm 1.5^\circ\text{C}$
Sulfate	$\leq 0.2\%$

Store at RT

PRODUCT USE LIMITATION
Research use only