

---

## PRODUCT SPECIFICATION

**Product Name:** RNase A, sol. 100 mg/mL

**CAT Number:** APL-0076

<b>Component</b>	
<b>Item</b>	<b>Concentration</b>
RNase A	0.073 mM

<b>Specification</b>	
<b>Format</b>	Ready to use solution
<b>Activity</b>	7,000 units/mL
<b>Unit definition</b>	One unit catalyzes the hydrolysis of RNA to yield a first-order velocity constant of 1.0 at 25°C, pH 5.0.

*Store at RT*

### PRODUCT USE LIMITATION

Research use only