


PRODUCT SPECIFICATION

Product Name: Tween 20

Chemical Structure:



CAT No.

APL-0086

CAS Registry No.

9005-64-5

Specification

DNase	None / Fail
RNase	None / Fail
Protease	None / Fail
Heavy Metals	<0.001%
Hydroxyl Number	96~108 mg KOH/g
Identification	Pass / Fail
Moisture (KF)	≤3.0%
Residue on Ignition	≤0.25%

Store at RT

PRODUCT USE LIMITATION

Research use only