


## PRODUCT SPECIFICATION

**Product Name:** TEMED

**Chemical Structure:**



**CAT Number:** APL-0092

**CAS Registry Number:** 110-18-9

**Molecular Formula:** C<sub>6</sub>H<sub>16</sub>N<sub>2</sub>

**Molecular Weight:** 116.20

<b>Specification</b>	
<b>Assay</b>	≥99.0%
<b>Appearance</b>	
<b>Color</b>	Colorless to Faint Yellow
<b>Form</b>	Liquid
<b>Clarity</b>	Clear
<b>Suitability</b>	Suitable for use in polyacrylamide gel electrophoresis

*Store at RT*

### PRODUCT USE LIMITATION

Research use only