
PRODUCT SPECIFICATION

Product Name: 10X Phosphate Buffer Saline (PBS), pH 7.3~7.5

CAT Number: APL-1041

Component	
Item	Concentration
KH ₂ PO ₄	14.70 mM
Na ₂ HPO ₄	95.50 mM
KCl	26.80 mM
NaCl	1.37 M

Specification	
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
pH at 25°C	7.3~7.5

Store at RT

PRODUCT USE LIMITATION
Research use only