

PRODUCT SPECIFICATION

Product Name: 20X Phosphate Buffer Saline Tween 20 (PBST) Wash Buffer, 1 L

CAT Number: APL-1043

Component	
Item	Concentration
Na ₂ HPO ₄ ·12H ₂ O	0.19 M
NaCl	2.74 M
KCl	0.05 M
KH ₂ PO ₄	0.03 M
Tween 20	1.00 %

Specification	
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
pH at 25°C	7.3~7.5

Store at RT

PRODUCT USE LIMITATION

Research use only