

AquaClean concentrate

CAT No. APL-5001

DESCRIPTION

AquaClean inhibits growth of bacteria, algae, fungi and precipitation of inorganic salts. AquaClean does not contain aldehydes or azides

CONTENT

	APL-5001
AquaClean concentrate	250 mL

STORAGE

Store at RT.

APPLICATION

Microbicidal additive for water soluble heating bath fluids.

PROTOCOL

1. AquaClean concentrate is supplied as a 500X concentrate and can be directly added to heating bath fluid at 1X (1:500). e.g. Add 2 mL AquaClean concentrate to 1 L of heating bath fluid.
2. Renew fluid and additive every 2 weeks.

PRODUCT USE LIMITATION

Research use only. Do not contaminate water, food, and feed by storage or disposal.